

**PATENT APPLICATION FOR UNITED STATES LETTERS PATENT  
IN THE UNITED STATES PATENT AND TRADEMARK OFFICE**

5

**(MBHB Case No. 02-1076-A)**

10    Title:            **Polymerized Hemoglobin Solutions Having Reduced Amounts of  
Tetramer and Method for Preparing**

15    Inventors:        Anthony Avella, a citizen of the United States and a resident of Crystal Lake, IL,  
Richard DeWoskin, a citizen of the United States and a resident of St. Charles, IL,  
Marc Doubleday, a citizen of the United States and a resident of Cary, IL.

20    Assignee:        Northfield Laboratories  
1560 Sherman, Suite 1000  
Evanston, IL 60201-4422

Express Mail No.: EV333549371US Date of Deposit: January 29, 2004
--

**Polymerized Hemoglobin Solutions Having Reduced  
Amount of Tetramer and Method for Preparing**

**CROSS REFERENCE TO RELATED APPLICATION**

[0001] This application claims the benefit of U.S. Provisional Application Serial. No. 60/443,436 filed January 29, 2003.

**FIELD OF THE INVENTION**

[0002] The invention relates to stabilized oxygen carrying solutions. More specifically, it relates to hemoglobin solutions that have been treated to enhance polymer bond stability and to remove tetramer that has been elaborated.

**BACKGROUND OF THE INVENTION**

[0003] There is a consistent need for ready blood products for an ever-increasing surgical and trauma load, and to supplement blood bank shortages. Oxygen carrying solutions, such as hemoglobin-derived solutions can be used in place of whole blood or red blood cells for patients having a need for augmented oxygen carrying capacity. Because they are not dependent upon donor availability, such solutions can be made readily available in an emergency situation or during a blood bank shortage. Also, due to risk of infection of blood borne pathogens as a result of a blood transfusion, a patient may prefer a hemoglobin-derived solution for transfusion in place of whole blood or red blood cells. In particular, such solutions may include, but are not limited to, oxygen carriers, blood substitutes, and hemoglobin-derived compositions such as

those disclosed in U.S Patent Nos. 6,498,141, 6,133,425, 5,464,814, 5,438,041, 5,217,648, 5,194,590, 5,061,688, and 4,826,811, the teachings of which are incorporated herein by reference in their entirety.

[0004] Stroma-free hemoglobin is known in the art to have oxygen transport and reversible  
5 oxygen (or ligand) binding capacities. However, hemoglobin solutions, while capable of carrying sufficient quantities of oxygen to support life, have presented challenges because of several undesirable side effects such as a decrease in kidney performance. These effects were thought to be due to the presence of unwanted contaminants such as bacterial endotoxin or fragments of red cell membranes (stroma) that is not removed from solution. While  
10 contaminants such as these can indeed produce renal alterations, hemoglobin solutions essentially free of such contaminants still produce substantial renal dysfunction. The cause for the renal dysfunction can be ascribed to, among other things, physiologically unacceptable amounts of unpolymerized hemoglobin tetramer.

[0005] Essentially tetramer free hemoglobin solutions can be used to replenish essentially all  
15 of a human patient's blood volume without causing vasoconstriction, renal toxicity, hemoglobinuria or other problems associated with intravenous administration of synthetic or semisynthetic oxygen carriers and blood substitutes. While these solutions provide superior efficacy, the shelf life of the product is limited since the hemoglobin polymer is known to slowly degrade to tetrameric units over time. Accordingly, what is needed is a method for maintaining

the solution as a substantially tetramer free solution for an extended period to increase the shelf life of the solution.

### **SUMMARY OF THE INVENTION**

[0006] In one aspect, the invention provides a method for producing a substantially tetramer  
5 free hemoglobin solution comprising polymerizing hemoglobin in solution, treating the  
polymerized hemoglobin to elaborate tetramer and removing tetramer from the polymerized  
hemoglobin solution.

[0007] In another aspect, the invention relates to a substantially tetramer free hemoglobin  
solution produced by polymerizing a solution of hemoglobin, heat treating the polymerized  
10 hemoglobin solution to elaborate tetramer, and removing tetramer from the polymerized  
hemoglobin solution.

[0008] In a further aspect, the invention relates to a method for stabilizing a polymerized  
hemoglobin solution comprising treating the polymerized hemoglobin solution to partially  
degrade the polymerized hemoglobin to tetramer and removing tetramer from the solution.

15 [0009] In yet another aspect, the invention provides a method for producing a stabilized,  
polymerized hemoglobin solution. The method includes producing a polymerized hemoglobin  
solution, removing tetramer from the polymerized hemoglobin solution to produce a  
substantially tetramer free polymerized hemoglobin solution, aging the polymerized hemoglobin  
solution to allow tetramer to elaborate, and removing the elaborated tetramer. The aging may

include storing the hemoglobin solution until the tetramer concentration is greater than about 1.0% of total hemoglobin, such as longer than one year.

[0010] The hemoglobin may be derived from mammalian blood, such as human or bovine blood. The hemoglobin may be polymerized with glutaraldehyde. The tetramer may be removed by filtration. The treatment of the polymerized solution to elaborate tetramer may be accomplished by heating the solution above about 45° C for at least about 24 hours. The tetramer concentration after removing tetramer may be less than about 1.0% of total hemoglobin in the solution, or less than about 0.3% of total hemoglobin in the solution.

10

#### **BRIEF DESCRIPTION OF THE FIGURES**

[0011] FIG. 1 shows the amount of tetramer elaboration in hemoglobin solutions stored at 2 – 8° C.

[0012] FIG. 2 shows the amount of tetramer elaboration in hemoglobin solutions stored at 23 – 27° C.

15 [0013] FIG. 3 shows a comparison of tetramer elaboration between newly processed hemoglobin solutions and reprocessed aged hemoglobin solutions stored at 2-8° C.

[0014] FIG. 4 shows a comparison of tetramer elaboration between newly processed hemoglobin solutions and reprocessed aged hemoglobin solutions stored at 23 – 27° C.

[0015] FIG. 5 shows a comparison of tetramer elaboration between newly processed hemoglobin solutions and hot quenched solutions stored at 2-8° C.

[0016] FIG. 6 shows a comparison of tetramer elaboration between newly processed hemoglobin solutions and hot quenched solutions stored at 23-27° C.

5 [0017] FIG. 7 shows a comparison of tetramer elaboration between newly processed hemoglobin solutions and pre-elaborated solutions stored at 2-8° C.

[0018] FIG. 8 shows a comparison of tetramer elaboration between newly processed hemoglobin solutions and pre-elaborated solutions stored at 23-27° C.

[0019] FIG. 9 is an HPLC tracing of a purified, polymerized hemoglobin solution.

10 Polymerized hemoglobin is indicated by the peaks at RT 15.19, 16.01, 17.16 and 18.79.  
Tetramer is indicated by peaks at RT 21.22 and 21.83.

[0020] FIG. 10 is an HPLC tracing of a polymerized hemoglobin solution after glycine treatment but prior to purification. Polymerized hemoglobin is indicated by peaks at retention times (RT) 14.62, 15.44, 16.60, and 18.24. Tetramer is indicated by the peak at RT 20.71.

15 Polymer is about 75 % of this material.

[0021] FIG. 11 is a schematic diagram depicting the portion of the process and equipment used to produce a deoxygenated hemoglobin solution prepared for pyridoxylation and polymerization.

[0022] FIG. 12 is a schematic diagram depicting the portion of the process and apparatus to produce, beginning with pyridoxylation and polymerization, a deoxygenated, purified, pyridoxylated, polymerized hemoglobin product, and the portion of the process and apparatus for formulating the final hemoglobin product having physiological levels of electrolytes.

5 [0023] FIG.13 is a schematic diagram depicting a column chromatography purification process employed in the invention.

[0024] FIG. 14 is a schematic diagram depicting a membrane filtration purification process employed in the invention.

10

#### **DETAILED DESCRIPTION**

[0025] The invention provides an oxygen carrying solution comprising an essentially tetramer-free, cross-linked, polymerized, hemoglobin which is substantially free of stroma, stromal contaminants and other contaminants.

[0001] Before describing the present invention in detail, a number of terms will be defined.

15 It is to be understood that the terminology used herein is for the purpose of describing particular embodiments only, and is not intended to be limiting. As used herein, the singular forms "a," "an," and "the" include plural referents unless the context clearly dictates otherwise.

[0026] "Hemoglobin" refers to hemoglobin derived from any source. Hemoglobin may be derived from mammals including humans, cattle, pigs, and sheep, or from other sources such as transgenically-produced hemoglobin described in BIO/TECHNOLOGY, 12: 55-59 (1994), and recombinantly produced hemoglobin, such as the recombinantly produced hemoglobin described  
5 in Nature, 356: 258-260 (1992). As used herein, % total hemoglobin (THb) is defined as grams of hemoglobin/100mL of solution.

[0027] A "solution of hemoglobin" refers to a solution of tetrameric or polymerized hemoglobin molecules where the molecules are not contained within a red blood cell. Such a solution need not be free of or substantially free of red blood cell stroma or stromal  
10 contaminants. However, in one aspect of the invention, solutions of polymerized hemoglobin are substantially free of red blood cell stroma and stromal contaminants.

[0028] "Cross-linked" means the chemical placement of molecular "bridges" onto or into a molecule, or between molecules with the purpose of altering the shape, size, function or physical characteristics of the molecule. Cross-linked molecules may be polymerized or non-  
15 polymerized, i.e., cross-linked molecules may be tetrameric.

[0029] "Tetramer" or "tetrameric" refers to hemoglobin molecules having a molecular weight of about 64Kd; that is, the term refers to both native and intramolecularly crosslinked hemoglobin molecules. As used herein % tetramer refers to the amount of tetramer as a



percentage of the amount of total hemoglobin (THb) in solution. For example, a 100 mL hemoglobin solution having 10% THb and 1% tetramer has 0.1 g tetramer in solution.

[0030] "Essentially tetramer free" denotes the level of purity with respect to tetramer contamination at which certain biological responses to tetramer administered into a mammal are no longer present. A main criterion is the absence of alterations in renal function when pharmaceutically effective amounts are infused, that is, at a level of purity of about 99% or better (less than about 1% of tetramer is present). The preferred product produced by the inventive process contains no more than about 0.8% tetramer based on the weight of total hemoglobin (THb). In other words, an essentially tetramer-free product according to the invention contains no more than physiologically acceptable amounts of unpolymerized hemoglobin tetramer. Particularly preferred products of the invention contain less than about 0.5% tetramer; the most particularly preferred products of the invention contain about 0.3-0.4% tetramer. Such amounts of tetramer have been found to be physiologically acceptable.

[0031] "Elaboration" or "tetramer elaboration" refers to an increase of the amount of tetramer in solution due to the degradation of polymerized hemoglobin to tetramer. Degradation of the polymer can be the result of treatment with chemicals, temperature, time or a combination thereof. Generally, tetramer elaborates upon heating or as the polymerized solution ages. Thus, percent tetramer increases during storage of the solution. Tetramer elaboration can be accelerated by heating the polymerized hemoglobin solution. Elaboration is said to have occurred after any increase in the amount of tetramer in solution. As used herein, the partial

degradation of a polymerized hemoglobin refers to some, but not all, of the polymer in solution being degraded to tetramer.

[0032] “Aging” a hemoglobin solution refers to storing the solution for any amount of time at any temperature. Higher temperatures accelerate the effects of aging on the hemoglobin solution. An “aged” hemoglobin solution has been stored so that tetramer has elaborated.

[0033] “Pre-elaboration” or “tetramer pre-elaboration” refers to a technique that utilizes heat treatment to promote tetramer elaboration.

[0034] “Hot quench” refers to a processing technique that involves heating the solution during the polymerization quench reaction to drive the reaction to completion.

[0035] “Polymerizing” or “Polymerized” means the act of, or the result of, the placement of molecular bridges between molecules or tetrameric subunits where the size and weight of the resulting polymerized molecule is increased with respect to native or tetrameric hemoglobin. Polymerized hemoglobin is not tetrameric hemoglobin. Polymerization may be accomplished using various polymerizing agents, including glutaraldehyde, imido esters, or others, in a biochemically suitable carrier, as is well known to those skilled in the art.

[0036] “Pyridoxylated” or “pyridoxylation” refers to the method of, or the result of, binding pyridoxal-5'-phosphate containing molecules to a hemoglobin molecule by reacting the molecule with pyridoxal-5'-phosphate (“P5P”) or 2-Nor-formyl pyridoxal-5'-phosphate. Pyridoxylation

has been shown to favorably alter the reversible oxygen binding capacity, i.e. increase the  $P_{50}$  of certain mammalian hemoglobins, *e.g.* human hemoglobin.

[0037] “Stable” or “stability” refer to the state or characteristic of hemoglobin solutions that are resistant to degradation and have a longer shelf life than non-stable solution. For example, 5 hemoglobin solution that have been stabilized according to the invention will, compared to solutions that have not been prepared according to the invention, have less or slower tetramer elaboration during storage of the solution. The stability of a hemoglobin solution is dependent on several other parameters that are independent of tetramer elaboration, including, for example, how quickly deoxyhemoglobin is converted to oxyhemoglobin or methemoglobin. This 10 parameter may be controlled by, among other ways, preventing oxygen from entering the packaged solution during storage. Stabilized hemoglobin solutions may still degrade, but do so at a slower rate than non-stabilized solutions.

[0038] The invention provides a polymerized, hemoglobin solution essentially free of tetrameric (native or intramolecularly crosslinked) hemoglobin, stromal and various other 15 contaminants. The solution is physiologically acceptable as well as therapeutically and clinically useful. The product has reversible oxygen binding capacity which is necessary for oxygen transport properties. The product demonstrates good oxygen loading and unloading characteristics in usage, which correlates to having an oxygen-hemoglobin dissociation curve ( $P_{50}$ ) similar to whole blood. The product binds oxygen with high affinity in the capillaries

through the lungs and then adequately releases oxygen to the tissues in the body. The product also does not require compatibility studies with the recipient.

[0039] In one aspect, the product may also have a half-life when administered to humans of about at least 15 hours. In another aspect, the half life is greater than about 24 hours. The  
5 hemoglobin product can be used to replenish essentially all of a human patient's blood volume without causing vasoconstriction, renal toxicity, hemoglobinuria or other problems associated with intravenous administration of synthetic or semisynthetic oxygen carriers and blood substitutes.

[0040] The half-life of the resulting product of the invention is determined *in vivo* in  
10 mammals, e.g., humans. Typically, a blood sample is removed from the mammal a period of time after the mammal has been infused with the product. The amount of the product is then determined by centrifuging the blood sample, expressing the plasma portion, determining plasma hemoglobin levels spectrophotometrically, and then correlating the amount of product remaining in the mammal to the half-life of the product.

15 [0041] The method of the invention yields a stabilized polymerized hemoglobin solution with an enhanced shelf life. Generally, tetramer elaborates as the polymerized hemoglobin solution ages. Elaborated tetramer can be removed from the aged hemoglobin solution by the process of the invention. Aged hemoglobin solutions that have been processed to remove tetramer exhibit slower tetramer elaboration upon further storage as compared to newly

processed (non-aged) hemoglobin solutions. Thus, hemoglobin solutions processed according to the invention exhibit a greater stability with respect to tetramer elaboration than newly processed solutions.

[0042] Generally, it has been found that physiologically and clinically useful hemoglobin solutions contain less than 1.0% tetramer. Accordingly, since tetramer elaborates over time during storage, to allow a reasonable shelf life for the hemoglobin solution, it is desirable that newly processed hemoglobin solutions contain less than about 0.3 % tetramer, with the expectation that tetramer in solution will elaborate, but the solution will remain physiologically useful until tetramer reaches 1.0%. FIG. 1 shows the typical amount of tetramer elaboration of a hemoglobin solution stored at 2 – 8° C. In a solution with beginning tetramer levels of about 0.5%, tetramer levels rise to greater than 1.0 % after about 18 months. FIG. 2 shows the typical amount of tetramer elaboration for a hemoglobin solution stored at 23 – 27° C. In a solution with a beginning tetramer level of about 0.5%, tetramer levels rise to above 1.0% in about 3-4 weeks. Thus, storage of hemoglobin solutions at increased temperatures increases the rate of tetramer elaboration. Aging of tetramer solutions can be accomplished merely by allowing the hemoglobin solution to stand for extended periods of time. Higher temperatures age the hemoglobin solutions faster than lower temperatures.

[0043] FIG. 3 shows a comparison of the tetramer elaboration between a newly processed hemoglobin solution and an “aged” hemoglobin solution processed according to the invention (“reprocessed”). Both solutions were stored at 2-8° C. After 24 months, tetramer increased

about 0.8% in the new hemoglobin solutions while, in the reprocessed solution, tetramer increased only about 0.3%. Over 36 months the tetramer in the reprocessed solution increased about 0.4%. Likewise, FIG. 4. shows that tetramer elaborates faster in newly processed solutions as compared to reprocessed solutions when stored at 23 – 27°.

5    **[0044]**    In one aspect, the method of the invention involves pre-elaborating the hemoglobin solution to enhances the stability of a solution by heating the solution to accelerate tetramer elaboration. The elaborated tetramer can be removed from the solution to provide a solution that is more stable with respect to tetramer elaboration than polymerized solutions that have not been heat treated to accelerate tetramer elaboration. Once elaborated tetramer has been removed from  
10   the hemoglobin solution, further tetramer elaboration is slowed or decreased.

**[0045]**    Heat treatment may occur either before or after purification of the solution. If the solution is heat treated during processing, such heat treatment should follow polymerization of the solution. The solution can then be purified to remove substantially all the tetramer. If heat treated after purification, the solution must be re-purified to remove the elaborated tetramer.

15   **[0046]**    Heat treating may be accomplished by subjecting the polymerized hemoglobin solution to about 45- 55° C for about 20-30 hours. It is expected that other processing temperatures and time will suffice in order to elaborate tetramer since it has been shown that tetramer elaboration is a function of, among other things, time and temperature.

[0047] The process of the invention provides a further advantage in that it can render the final product substantially free of microbial and viral antigens and pathogens. Such microbial and viral antigens and pathogens are reduced to non-detectable levels, i.e., the product is sterile as determined by the analysis set forth in the United States Pharmacopoeia, XXIII Chapter 71.

5 Examples of such antigens and pathogens include, for example, bacterial, rickettsial, fungal, protozoan, viral and other organisms. Most importantly, the process provides a biological product free of viruses that cause hepatitis and acquired immune deficiency syndrome (AIDS).

[0048] Heat treatment of the hemoglobin solution results in the substantial inactivation of viruses. Viral inactivation heat treatment can be a separate step of heat treatment or may be  
10 combined with the heat treatment to remove tetramer. Generally, for the purposes of production operator safety, it is preferred that viral reduction heat treatment be conducted after the hemoglobin has been removed from red blood cells. However, it is expected that heat treatment of the solution following polymerization will also accomplish the desired viral inactivation. If a pre-polymerization viral reduction heat treatment step has been eliminated, temperature of the  
15 post-polymerization heat treatment to elaborate tetramer can be increased to about 60-62° C, or other suitable temperature, to ensure viral inactivation. Thus, one heat treating step following polymerization will, in certain aspects, accomplish both viral reduction and tetramer elaboration.

[0049] The biological product of this invention, when infused in amounts of up to at least about 10.0 L, does not cause vasoconstriction, renal toxicity, hemoglobinuria and other problems  
20 implicated with intravenous administration of known hemoglobin solutions containing

physiologically undesirable amounts of tetrameric hemoglobin. Intravenous administration of the product produced by the process described herein results in no appreciable decrease in urine production, no appreciable decrease in glomerular filtration rate, no appreciable extravasation into the peritoneal cavity and no appreciable change in the color of urine produced.

5   **[0050]**     Therefore, the process of the invention provides an acellular red blood cell substitute useful in the treatment of trauma, myocardial infarction, stroke, acute anemia and oxygen deficiency disorders such as hypoxemia, hypoxia or end stage hypoxia due to impairment or failure of the lung to fully oxygenate blood. The product also is useful in the treatment of any disease or medical condition requiring a resuscitative fluid (*e.g.*, trauma, specifically  
10 hemorrhagic shock), intravascular volume expander or exchange transfusion. In addition to medical treatment, the product can be useful in preserving organs for transplants.

**[0051]**     In one aspect, the starting material in the process of the invention is whole human blood or packed red blood cells. Generally, it is desirable, but not critical, to use source red blood cells that have been in storage for no more than 2 weeks past the expiration date indicated  
15 on the blood storage bag. The use of human whole blood outdated by more than 2 weeks provides additional difficulty in extracting the hemoglobin and removing cellular remnants such as stromal proteins and contaminants. In addition, the processes described herein are applicable to all hemoglobins with minor modifications within the skill of the art.



[0052] If human blood is used as a starting material, during red cell aspiration and filtration, the red blood cells (RBC) are aseptically extracted from donor bags without introducing air into the blood and passed across a series of filters to result in a RBC suspension having reduced amounts of leukocytes and platelets. The resulting suspension is then subjected to cell  
5 washing/lysing.

[0053] The suspension is washed under carbon monoxide atmosphere with an about 1% NaCl solution to remove residual plasma proteins. The washed RBC are then treated with water for injection ("WFI") to lyse the cells and the resulting mixture is clarified using a cross flow filtration unit. Other methods of lysing red blood cells known to those of skill in the art may be  
10 used including, for example, mechanically or sonically lysing the cells. The clarified product may then be heat-treated for viral inactivation and to precipitate additional stromal material which is removed by filtration. The product of this procedure is a stroma-free hemoglobin (SFH) solution with a THb of about 3% (w/v).

[0054] Following clarification, the solution containing carboxyhemoglobin is preferably  
15 concentrated and degassed to yield a stroma free hemoglobin solution containing deoxyhemoglobin. Degassification involves first saturating the carboxyhemoglobin solution with oxygen for about 16 hours to yield a solution of oxygenated hemoglobin and about 7% by weight, based on the total weight of hemoglobin, of carboxyhemoglobin. Subsequently, the oxygen is driven off with nitrogen, argon or helium to form a solution containing free  
20 hemoglobin, i.e., uncomplexed hemoglobin, and about 7% by weight, based on the total weight

of hemoglobin, of oxyhemoglobin. The resulting degassed solution is filtered and transferred into a vessel for chemical modification.

[0055] Subsequent to degassification, the stroma-free hemoglobin solution comprising human hemoglobin should be pyridoxylated using pyridoxal-5'-phosphate (P5P) at a molar ratio of pyridoxal-5'-phosphate to hemoglobin of about 1:1 to 3:1. Alternatively, the stroma-free hemoglobin may be pyridoxylated using 2-Nor-2 formyl pyridoxal-5'-phosphate. A reducing agent such as sodium cyanoborohydride or preferably sodium borohydride is added to the pyridoxylation mixture. Excess reagents and salts are removed by dialysis against pyrogen free water or, preferably, diafiltration with WFI. Hemoglobin from sources other than human blood may not need pyridoxylation. Those skilled in the art of hemoglobin solutions readily understand when pyridoxylation is required.

[0056] The stroma-free, hemoglobin solution is polymerized using any method known to those skilled in the art of hemoglobin solutions. Preferably, an aqueous glutaraldehyde is used as a polymerizing agent. The duration of polymerization and the amount of glutaraldehyde added is dependent on volume of the hemoglobin solution, the desired yield of polymers and the desired molecular weight distribution. In general, longer polymerization times increase the yield and the molecular weight distribution of the polymers. A yield of approximately 75% by weight of polymers, based on the total weight of hemoglobin, is obtained in about 16-18 hours. The preferred end point of the polymerization is defined as that point where the solution contains about 75% by weight of polymers, based on the total hemoglobin weight, as monitored by size-

exclusion HPLC. Alternatively, the endpoint is defined as the point at which the solution contains about 65% of polymers based on the total weight of hemoglobin, i.e., about 2.5 hours.

[0057] Following polymerization, the reaction should be quenched with the appropriate agent. In one aspect, the polymerization reaction is quenched by the addition of aqueous glycine.

5 The glycine should be added as quickly as possible. The cross-links are then stabilized by adding, again as quickly as possible, a solution of aqueous sodium borohydride. This polymerized solution is subsequently concentrated and then diafiltered. Water is finally added to the solution until the solution contains about 4% by weight hemoglobin.

[0058] In another aspect, the solution may be "hot quenched" by heating the solution to 40-  
10 50° C for at least three hours concurrent with the addition of glycine to drive the quench reaction to completion. FIG. 5 shows a comparison of the tetramer elaboration between a newly processed hemoglobin solution and a hemoglobin solution subject to a hot quench for three hours. Both solutions were stored at 2-8° C. After 10 months, tetramer increased between 0.4-0.5% in the new hemoglobin solution while, in the hot quenched solution, tetramer increased  
15 only about 0.4%. FIG. 6 shows the differences in tetramer elaboration between newly processed solutions and hot quenched solutions when stored at 23 – 27° C. It is expected stability of the solution would be enhanced by subjecting to the solution to a longer hot quench, for example, up to twenty-four hours. Higher temperatures, up to about 65° C, may also be employed.

[0059] Polymerization according to the invention results in a high yield of polymers having a narrow molecular weight range as shown in FIG. 9 and the Examples below.

[0060] In another aspect, the post-quench polymerized solution may be pre-elaborated by heat treatment to elaborate tetramer. The heat treatment may be postponed until any point after purification, but the solution would require repurification. Heat treatment may be accomplished by heating the solution to above about 45°C for at least about 24 hours. If viral inactivation is desired at this point, the solution may be heated to above about 60°C. An antioxidant such as ascorbic acid may be added to prevent formation of methemoglobin during the heat treating process.

[0061] FIG. 7 shows a comparison of the tetramer elaboration between a newly processed hemoglobin solution and a pre-elaborated hemoglobin solution processed according to the invention. Both solutions were stored at 2-8° C. After 15 months, tetramer increased about 0.5% in the new hemoglobin solutions while, in the pre-elaborated solution, tetramer increased only about 0.4%. Likewise, FIG. 8 shows that tetramer elaborates faster in newly processed solutions as compared to pre-elaborated solutions when stored at 23 – 27° C.

[0062] The polymerized, pyridoxylated hemoglobin solution is then purified. In one aspect, purification is accomplished under an atmosphere of oxygen to oxygenate the solution utilizing column chromatography, filtration, e.g., membrane filtration, or both, to remove residual unpolymerized (tetrameric) hemoglobin from the solution. The purified polymerized

hemoglobin solution is then concentrated to about 6% using an ultrafiltration apparatus in preparation for gas exchange.

[0063] The concentrated solution is then deoxygenated with nitrogen. The deoxygenation takes place at about 10-12°C until the amount of oxyhemoglobin in the solution is less than about  
5 16% by weight of the total hemoglobin.

[0064] The resulting deoxygenated, purified, and polymerized hemoglobin solution is then concentrated by ultrafiltration under a nitrogen atmosphere in a cooled vessel. The pH is adjusted to about 8.8-9.0, and the amounts of electrolytes may be adjusted as necessary to levels representing that of normal plasma. In addition, conventional antioxidants such as glutathione,  
10 ascorbate or glucose may be optionally added. After the solution is concentrated to the desired level, preferably about 10% by weight total hemoglobin, the solution is sterilized by filtration and transferred via a sterile transfer apparatus into suitable pharmaceutically acceptable containers.

[0065] If the hemoglobin solution was not previously heat treated to facilitate tetramer  
15 elaboration, the solution may be heat treated and repurified to remove elaborated tetramer. If antioxidants and formulation chemicals have been added, these may be removed by diafiltration prior to heat treatment and purification.

[0066] In addition, as an alternative to heat treatment, tetramer elaboration of the finished hemoglobin solution may be accomplished by allowing the solution to age at an appropriate

temperature. The solution can then be repurified to remove the elaborated tetramer. Generally, it is expected that the solution should age until tetramer levels exceed about 1 – 3 %, with higher tetramer levels giving increased stability benefits after the solution is purified to remove tetramer. However, it is expected that the advantages of the invention will be accomplished if the solution is aged for any period of time, so long as the elaborated tetramer is removed after aging of the solution.

[0067] The characteristics of the resulting hemoglobin solution are shown in FIG. 9 and are as follows:

Total Hemoglobin (g/dl)	9.5-10.5
Methemoglobin (% of total Hb)	< 8.0
Carboxyhemoglobin (% of total Hb)	< 5.0
P <sub>50</sub> (torr)	26-32
Osmolality (mmol/Kg)	280-360
Sodium (mmol/L)	135-155
Potassium (mmol/L)	3.5-4.5
Chloride (mmol/L)	85-110
Free Iron (ppm)	< 2.0
Molecular Wt. Dist. – 128 Kd peak (%)	9-23
Molecular Wt. Dist. – 192 Kd peak (%)	16-18
Molecular Wt. Dist. – 256 Kd peak (%)	49-74
Tetramer (64K)(%)	≤ 1.0
Endotoxin (EU/mL)	< 0.03
Phospholipids ng/Hb	< 50

Glycolipids (ng/Hb)

< 2

### Examples

[0068] The following examples demonstrate certain aspects of the invention. However, it is to be understood that these examples are for illustrative purposes only and do not purport to be wholly definitive as to conditions and scope of this invention. All temperatures are expressed in degrees Celsius unless otherwise specified. It also should be appreciated that when typical reaction conditions (e.g., temperature, reaction times) have been given, the conditions which are both above and below these specific ranges can also be used, though generally less conveniently.

[0069] Unless noted to the contrary, all vessels and tanks used in the inventive process are made of 316-L Stainless Steel, preferably a pharmaceutical grade of such stainless steel that has been highly polished and therefore easily and rapidly cleaned. The various connecting pipes and tubes are made of the same stainless steel or of a pharmaceutical grade Teflon or silicone tubing. The filters and membranes used in the process may be purchased from Millipore Inc., Pall-Filtron, or Cuno Inc.

[0070] Analytical Size Exclusion Chromatography HPLC according to the invention is carried out according to the following procedure. The sample is diluted to 0.2 g/dl with 0.1 M sodium phosphate buffer at about pH 6.9, filtered through a 0.2 $\mu$  filter and injected into an HPLC system consisting of the following components (in order of system flow):

1. Agilent Technologies 1100 Isocratic Pump

- mobile phase is 0.1 M sodium phosphate at about pH 6.9
- flow rate is 0.5 mL/minute
- 2. 45 cm PEEK or titanium tubing, 0.010 in. I.D.
- 3. Agilent 1100 Autosampler
- 5 4. 18 cm PEEK or titanium tubing, 0.010 in. I.D.
- 5. 0.5  $\mu$  precolumn filter frit
- 6. 9 cm PEEK or titanium tubing, 0.010 in. I.D.
- 7. Toso TSK G3000SWXL 40 x 60 mm guard column
- 8. 24 cm PEEK or titanium tubing, 0.010 in. I.D.
- 10 9. Toso TSK G3000SWXL 300 x 7.8 mm Analytical column
- 10. 23 cm PEEK or titanium tubing, 0.010 in. I.D.
- 11. Agilent 100 variable wavelength detector.
  - wavelength: 280 nm
  - flow cell: 14  $\mu$ L vol., 10 mm pathlength
  - 15 - range: 2 AUFS
  - time constant: 10 seconds

The peak absorbance at 280 nm is recorded by an Agilent Chemstation, which integrates the individual peak areas and calculates the total Hemoglobin area for each polymeric species.



## **EXAMPLE 1**

### **Cell Aspiration and Filtration**

[0071] Referring now to FIG. 11, donor bags 20 of outdated blood (whole blood or packed red blood cells) are situated in a suitable aseptic aspiration apparatus 22. A needle in the aspiration apparatus punctures the donor bag, introduces about 150 ml of a 1% (w/v) aqueous sodium chloride solution and aspirates the outdated blood from the donor bag under reduced pressure or vacuum. The aspirated blood is passed through leukocyte adsorption depth filter 24 or alternatively through two 5 $\mu$  depth filters in series 26. As the blood passes through the filters, leukocytes are removed from the blood. Typically, about 225 units of outdated whole blood are aspirated, filtered and subsequently transferred to Tank 1 as shown in FIG. 11. The filters are then rinsed with about 75 liters of a 1% (w/v) aqueous sodium chloride solution.

## **Example 2**

### **Cell Wash and Lysis**

[0072] Prior to the introduction of the blood into Tank 1, Tank 1 is charged with about 40-50 L of a 1% aqueous sodium chloride solution. After all 225 units of outdated whole blood have been aspirated, filtered and transferred, and the filters have been rinsed, the tank contains about 365-395 liters of a 4% total hemoglobin solution. During the aspiration and filtering steps, Tank 1 is maintained at a reduced pressure, i.e., a vacuum of 20-28 inches Hg. Once all the outdated blood has been transferred to Tank 1, the vacuum is switched off and carbon monoxide is introduced into the tank so that the tank contains an atmosphere of carbon monoxide.

[0073] Tank 1 is coupled to a 0.65μ tangential flow filter 28 as shown in FIG. 11. The initial charge of 365-395 liters of 4% total hemoglobin solution is concentrated to approximately 215-225 L of a 7% total hemoglobin solution by microfiltration through the tangential flow filter.

The pH of the hemoglobin solution at this point is about 6 to 6.5. Subsequent to concentrating to

5 7% total hemoglobin, the solution is washed by adding a 1% (w/v) sodium chloride solution, diafiltering and removing the filtrate at the same rate sodium chloride solution is added. The 215-225 L of hemoglobin solution is typically washed with about 8 volumes of the 1% sodium chloride solution (about 1,800 L). Subsequent to washing, the solution is concentrated to about 90-95 L, i.e., about 16% total hemoglobin, and water for injection ("WFI") is added to bring the  
10 volume of the solution up to about 220 L. With the addition of the WFI, the cells swell and rupture releasing hemoglobin into solution. The concentration of the resulting hemoglobin solution is about 7% total hemoglobin (THb).

[0074] The resulting solution is clarified while still in Tank 1. The solution is first

concentrated to about 90 L and the filtrate is transferred to Tank 2. As the solution is pumped

15 across the filter, red blood cells, stroma contaminants and cell wall material is retained and removed by the filter. The remaining 90 L of solution in Tank 1 is washed (diafiltered) with about 280 L of WFI and the wash is added to Tank 2. The material remaining in Tank 1 is then concentrated to about 20 L and the filtrate added to Tank 2. The volume resulting in Tank 2 is about 405-415 L of a 3.3% total hemoglobin solution.

### **Example 3**

#### **Optional Heat Treatment For Viral Reduction and Stromal Precipitation**

[0075] The resulting solution of stroma-free hemoglobin is then heat treated in Tank 2 at a temperature of about 60-62°C over a period of about 10 hours. During this time, the solution is moderately agitated. As the solution is heated and passes a temperature of about 55°C, a precipitate forms.

### **Example 4**

#### **Clarification and Viral Reduction**

[0076] The resulting 3.3% THb stroma-free, heat treated hemoglobin solution is then filtered through a 0.2μ prefilter 30 followed by a 0.1 μ prefilter 32 and then pumped through a 100 kD viral reduction ultrafilter (not shown) into Tank 3.

### **Example 5**

#### **Ultrafiltration/Concentration**

[0077] The filtered hemoglobin solution is then concentrated to 85-105 L (about 14% THb) and subsequently washed and diafiltered with about 4 volumes of WFI (350 L). The concentration and diafiltration is accomplished using a 10 kD molecular weight ultrafilter 34. Drain 35 associated with ultrafilter 34 collects filtrate. At this point, the 14% total hemoglobin solution contains less than about 50 ng of phospholipid per gram of hemoglobin, less than about 2 ng of glycolipid per gram of hemoglobin, less than about 1% methemoglobin, less than about

0.03 units of endotoxin per milliliter at a pH of about 6 to 6.5. The hemoglobin in the solution is carboxyhemoglobin.

## **Example 6**

### **Degassification**

5   **[0078]**   The resulting carboxyhemoglobin solution is then transferred to a 175L vessel (Tank 4) where the carboxyhemoglobin is first oxygenated and then deoxygenated. Tank 4 is fitted with a gas sparge ring coupled to oxygen and nitrogen gas lines, a feed from the tank bottom to a metered spray apparatus positioned at the top of Tank 4, and a foam overflow collector connected to foam can 36 such that foam generated in Tank 4 is fed into foam can 36 where the  
10   foam condenses into liquid and is fed back into Tank 4. As an alternative to the foam can 36, Tank 4 can be fitted with a mechanical foam breaker. Tank 4 further includes a center mounted, gas dispersion agitator. Foam can 36 includes a gas vent for removal of gas. The solution in Tank 4 is a 13% total hemoglobin solution.

**[0079]**   During a first oxygenation step, oxygen is sparged through the solution at a rate  
15   sufficient to have uniform dispersion of gas in the vessel. The vessel is sparged at a rate of about 66 L/min. with gas. Oxygenation of the carboxyhemoglobin is conducted for a period of about 18 hours such that the resulting solution contains less than 5% carboxyhemoglobin based on the weight of total hemoglobin. Oxygenation is conducted at a temperature of about 10°C. The foam generated in Tank 4 is collected in Foam Can 36 and after settling, the resulting solution is  
20   transferred back into Tank 4.

[0080] After oxygenation, the solution is sparged with a similar flow of nitrogen for about 3-3.5 hours or until less than 10% oxyhemoglobin based on the weight of total hemoglobin remains in the solution. The nitrogen sparge is conducted at a temperature of about 10°C and a pH of about 6.95 – 7.10. Alternatively, carboxyhemoglobin could be converted to deoxyhemoglobin using a membrane exchanger. It is noted that there is substantially no denaturing of the hemoglobin as would normally be expected from the foaming step.

### **Example 7**

#### **Chemical Modification**

[0081] Referring now to FIG. 12, the deoxyhemoglobin solution is transferred to Tank 5 for chemical modification. To Tank 5 containing the deoxyhemoglobin solution at about 4°C is then added an aqueous solution of pyridoxyl-5-phosphate (P5P) (93.75 g/L) at a 1:1 to 3:1 P5P to hemoglobin molar ratio. A 2:1 molar ratio of P5P to hemoglobin is preferred. The pyridoxylation is conducted at a temperature of about 4°C. The P5P solution is typically added over about 1 minute and mixed for approximately 15 minutes, after which a sodium borohydride/sodium hydroxide solution is added to the hemoglobin solution at a molar ratio of sodium borohydride to hemoglobin of about 20:1. A suitable aqueous sodium borohydride/sodium hydroxide solution contains 0.8 g of sodium hydroxide per 2 liters and 90.8 g of sodium borohydride per 2 liters. The borohydride solution is added as rapidly as possible over a period of about 1 minute and then stirred for one hour.

## **Example 8**

### **Reactant Removal**

[0082] The resulting 150 L solution of pyridoxylated hemoglobin is subsequently diafiltered using 10K Dalton ultrafilter 38 to remove excess reactants with 4 volumes (600 L) of WFI.

- 5 Drain 40 associated with ultrafilter 38 collects the filtrate from filter 38.

## **Example 9**

### **Polymerization**

[0083] To Tank 5 containing the pyridoxylated hemoglobin is added sufficient WFI to prepare a 4.5% total hemoglobin solution (about 265 L of hemoglobin solution). A

- 10 glutaraldehyde solution is added to the pyridoxylated hemoglobin solution at a molar ratio of glutaraldehyde to hemoglobin of about 24:1. The glutaraldehyde solution is typically added over a period of about 2.5 hours by a metering pump to the hemoglobin solution. The polymerization reaction is allowed to proceed for about 18 hours. The target molecular weight distribution is about 75% polymer and 25% tetramer. The target polymers have molecular weights of less than  
15 about 600,000 with a predominant fraction of the molecular weights residing in the 100,000-350,000 range.

[0084] When the polymerization reaction reaches the target molecular weight distribution (after about 18 hours), aqueous glycine (about 166 g/L) is added (as a quench) to the hemoglobin solution at a 140:1 molar ratio of glycine to hemoglobin. At this point, the solution may be

heated to 40-50° C for at least three hours to drive the quench reaction to completion (“hot quench”).

[0085] FIG. 10 shows an HPLC tracing of the resulting polymerized, glycine-quenched hemoglobin product. The resulting solution is then mixed for about 10 minutes after which a sodium borohydride sodium/hydroxide solution (having the concentration identified above) is added to the hemoglobin solution at a 28:1 molar ratio of sodium borohydride to hemoglobin. This resulting mixture is stirred for about 1 hour. The solution is then concentrated to about 150 L (ultrafilter 38) and washed with 4 volumes (600 L) of WFI. An additional aliquot of sodium borohydride at the same molar ratio as indicated above is added to the concentrated solution and again mixed for 1 hour. The resulting solution is washed with 4 volumes of WFI (600 L).

#### **Example 10**

##### **Optional Heat Treatment for Tetramer Elaboration**

[0086] At this point, the solution may be heat treated to elaborate tetramer. The solution may be subject to heating at 45 – 55 °C for about 20 – 30 hours. If the viral reduction is desired during this step, the temperature may be increased to above 60° C. To heat the solution, a heating medium, such as a propylene glycol solution at about 80°C, is circulated through the tank jacket while vigorously agitating the Hb solution. After heating, the hemoglobin solution is cooled to about 2-8°C.

## Example 11

### Purification

[0087] The resulting solution is oxygenated by allowing the solution to stand under an oxygen atmosphere and is subsequently transferred to a purification system 42. The purification  
5 may be achieved by column chromatography, filtration, preferably membrane filtration (diafiltration), or a combination of filtration and column chromatography.

[0088] In one embodiment, the solution is transferred to chromatography feed vessel, Tank 6, as shown in FIG. 13. In this embodiment, the resulting solution of oxyhemoglobin is then diluted to about 200 L (4% total hemoglobin) in Tank 6 and the concentration of chloride is  
10 adjusted to 22 mM with sodium chloride solution. No adjustment of sodium concentration is necessary.

[0089] Five 40 L aliquots of the resulting hemoglobin solution are then chromatographed using Column 44. Column 44 contains an affinity gel which is an agarose gel modified with a yellow dye (commercially available from Affinity Chromatography, Ltd., as Mimetic Yellow  
15 No. 1) having greater affinity for polymer than tetramer.

[0090] The chromatography is accomplished as follows. 40 L of oxygenated, polymerized, pyridoxylated, stroma-free hemoglobin solution is loaded onto Column 44. The column is washed with 15 column volumes (about 750 L) of 30 mM aqueous NaCl buffer to remove tetramer. The column is then washed with about 250 L of a 300 mM sodium chloride buffer to  
20 wash the polymer off. Polymer fractions are collected in Tank 7. Unwanted fractions are sent



to drain 46. After each aliquot is removed, the column is regenerated with 15 mM HCL solution (150 L), re-equilibrated with 30mM aqueous NaCl (250 L) and another aliquot of feed solution (40 L) is loaded to the column. The column is again washed with 30 mM NaCl followed by 300 mM NaCl. 40 L aliquots of hemoglobin solution are added to the column and chromatographed  
5 until Tank 6 is empty.

[0091] The collected fractions in Tank 7 are ultrafiltered (concentrated) using filter 48 associated with drain 50 to a volume of about 40 L (6% total hemoglobin). The concentrated hemoglobin solution is then transferred to gas exchange Tank 8 for deoxygenation.

[0092] Alternatively, the solution is transferred to a filtration recycle vessel 10, as shown in  
10 FIG. 14. The hemoglobin is then diluted to about 4% THb in Tank 10. The 4% THb solution is then diafiltered using 10 mM NaCl and a 300,000 molecular weight filter 52 commercially available from Millipore Corporation. The filtration is continued until about 97% of the hemoglobin material passes through the filter and into Tank 11. (About 3% of the material, i.e., high molecular weight polymers, is retained in Tank 10). The amount of hemoglobin is  
15 determined spectrophotometrically using a cooximeter.

[0093] The resulting material in Tank 11 is about 4-8% THb and contains about 7-10% tetramer based on THb. The 4-8% THb is then diafiltered using 10mM NaCl and a 100,000 molecular weight filter 54 commercially available from Pall-Filtron associated with drain or trap  
56. The filtration is continued until the level of tetramer, as determined by size exclusion

chromatography using a Toso BioSep 300 x 7.8 mm column, is less than 1.0% of the hemoglobin mass by weight. The resulting purified hemoglobin solution remains initially in Tank 11 and is subsequently transferred to gas exchange Tank 8 for deoxygenation.

## **Example 12**

### **5 Deoxygenation**

[0094] Gas exchange Tank 8 may be the same tank as Tank 4 or, preferably, a different tank. Gas exchange Tank 8 is equipped in essentially the same fashion as gas exchange Tank 4 and is attached to foam can 58 or equipped with a mechanical foam breaker in a fashion identical to that of Tank 4. Deoxygenation is accomplished in about 2.5 hours with a nitrogen sparge at  
10 about 10°C and a solution pH of about 8.8. Nitrogen sparging is continued until less than about 16% oxyhemoglobin, based on the weight of total hemoglobin, remains in the solution. The resulting deoxyhemoglobin solution is subsequently transferred to Tank 9 for formulation.

[0095] If the heat treatment step of Example 10 has been postponed, it may be conducted at this stage. If heat treatment is conducted now, the purification procedures of Examples 11 and  
15 the deoxygenation procedures of this Example should be repeated. The heat treatment step of Example 10 and the heat treatment step at this point are optional and may be eliminated completely.

### Example 13

#### Formulation

[0096] In formulation Tank 9, the solution is first concentrated to about 7% total hemoglobin, and the pH is adjusted to about 8.8 to 9.0 at 4°C. The pH is adjusted using 0.2 M NaOH. Glucose and glycine are added to achieve final concentrations of about 1 g/L and 3.5 g/L respectively. Potassium chloride is added to the solution to obtain a potassium concentration of about 3.5 to 4.5 mM. 3 M sodium chloride is then added to obtain an 85-110 mM chloride concentration. Sodium lactate is subsequently added to obtain a 135-155 mM concentration of sodium ion. Finally, a 0.45 molar ascorbic acid solution is added until the ascorbic acid concentration reaches about 1000 mg/L. The pH is adjusted to 8.7-9.1 at 10-15° C using 0.22 M NaOH. The resulting hemoglobin solution has a final osmolality of about 280-360 mmole per kg.

[0097] The formulated hemoglobin solution is then concentrated to about 10% total hemoglobin using filter 60 associated with trap 62. The 10% hemoglobin solution is then sterilized by filtration through filter 64 and aseptically filled into presterilized bags.

### Example 14

#### Solution Pre-Elaboration

[0098] Whether or not the solution has been subjected to the hot quench or optional heat treatment steps to elaborate tetramer during processing, the solution may be heated anytime after formulation to elaborate tetramer by heating the solution to about 45 – 55° C for about 20-30

hours, or until tetramer has increased above about 1 - 3 %. Following heat treatment the solution should be purified as in Example 11 and degassed as in Example 13. Since the solution at this point has already been formulated, the additives introduced in Example 12 are preferably removed by diafiltration prior to purification, and then reintroduced during re-formulation.

5 **Example 15**

**Re-processing Aged Hemoglobin Solutions**

[0099] Out of date or aged hemoglobin solutions may be reprocessed as in Example 14, without heat treating. Typically, hemoglobin solutions that are older than about 12 – 18 months at storage conditions of about 2-8°C have tetramer levels of greater than about 1.0%. This material may be introduced to Tank 5 and purified and reformulated according to Examples 11 - 13. Since the solution at this point has already been formulated, the additives introduced in Example 13 are preferably removed by diafiltration prior to purification, and then reintroduced during re-formulation.

[00100] In the foregoing, there has been provided a detailed description of preferred embodiments of the invention for the purpose of illustration and not limitation. It is to be understood that all other modifications, ramifications and equivalents obvious to those having skill in the art based on this disclosure are intended to be within the scope of the invention as claimed.